




## TECHNICAL SPECIFICATIONS

Type	Optical imaging platform for ocular diagnosis
Image resolution	Up to 3088 x 2076 px
Acquisition mode	Photo or video
Light source	White LED, Blue LED, Infrared LED
Yellow filter	Yes
Focus	Automatic or manual
Exams	Tear stability, Blinks, Lipid layer, Tear meniscus, Standard Meibography and Trans-illumination Meibography, Fluo stability, Generic exams (in colors and in fluorescein)
Questionnaires	Triage (TFOS-DEWS II), symptoms (OSDI©, SPEED™, DEQ-5)
Grading scales	Efron, Oxford, Guillon, Pult
DICOM	Import patient worklist
Connection	4 USB ports, 1 Ethernet port, 1 HDMI port, Wifi, Bluetooth
Screen	Color touch screen 15"
Computer	1 To SSD hard disk, Windows® 10 Enterprise
Eyelid everter	Standard Meibography and Trans-illumination Meibography
Temperature	15-30 °C
Power supply	110/230 VAC, 50/60 Hz, 300 VA
Dimensions	Standard Configuration: 500 (W) x 500 (H) x 500 (D) mm Configuration with secondary screen: 600 (W) x 500 (H) x 500 (D) mm
Weight	15 kg
Accessories in option	Secondary screen for exam visualization by the patient (10.1"), examination table with keyboard drawer

Specifications are subject to change without notice. Non contractual pictures.  
© 2023, C.Diag® is a trademark of Quantel Medical and Lumibird Medical. All rights reserved.  
XS\_CDIAG\_BC1\_AN\_230818

### 3 CONFIGURATIONS AVAILABLE

#### ① Configuration with dedicated examination table:



#### ② Standard configuration:



#### ③ Configuration with screen for exam visualization by the patient:



**Product not FDA approved & not CE approved**

#### Manufacturer

Quantel Medical  
1 Rue du Bois Joli - CS40015  
63808 Cournon d'Auvergne - FRANCE  
Tel. : +33 (0)4 73 745 745  
Email: [contact@lumibirdmedical.com](mailto:contact@lumibirdmedical.com)  
ISO 13485 : 2016

#### Headquarters

Lumibird Medical  
1 Rue du Bois Joli - CS40015  
63808 Cournon d'Auvergne - FRANCE  
Tel. : +33 (0)4 73 745 745

**QUANTEL MEDICAL**  
A brand of the group



**LUMIBIRD®**  
MEDICAL

[www.lumibirdmedical.com](http://www.lumibirdmedical.com)